


**New!**

# The Spectrum<sup>CS</sup> In-Clinic System

## A Breakthrough for Routine Culture and Sensitivity Testing

- Easily interpret color reactions for Staphylococcus, E. coli, Pseudomonas, and other common organisms.
- Differentiate Gram (+) from Gram (-) organisms without staining.
- Incorporates gold-standard disc diffusion method for antibiotic susceptibility testing.
- Comprehensive results in as little as 18-24 hours.
- Economical and efficient.
- Applicable for all species.



Spectrum™ CS is the complete in-clinic solution for identifying and determining antimicrobial sensitivity for many common bacterial organisms. At the heart of the Spectrum CS system is the new Spectrum IV quad plate that combines Gram positive and Gram negative chromogenic agars, selective Staphylococcus agar, and traditional blood agar into one convenient product.

What makes Spectrum™ chromogenic agars so innovative is that each has been formulated with special substrates that produce colored colonies unique to each bacterium. As an example, E. coli organisms produce pink to red colonies while Klebsiella appear blue to dark blue. Pseudomonas colonies are light green. Streptococci appear as pinpoint light blue colonies.

Incorporated into the Spectrum CS system is a separate Mueller Hinton plate, the medium of choice for antibiotic susceptibility testing. Also included is the BBL Prompt system that simplifies preparation of standardized organism suspensions. Antibiotic sensitivity discs (sold separately) can be selected from a menu of over 40 antimicrobials allowing you to choose only those that fit the specific needs of your practice.




*BBL Prompt® System*

The Spectrum CS kit includes enough plates and accessories to process 5 specimens. All plates come individually wrapped for extended shelf life.

*Distributed exclusively by:*



18131 SW 98th Ct.  
Palmetto Bay, FL 33157  
800.330.1522 · 305.232.8421 Fax  
www.vetlab.com



**Full color laminated interpretation guide provided free with first order.**



**New and refurbished incubators tailored to your test volume and budget.**

**Ordering Information:**

- PLTSP500      Spectrum CS Kit; 5/Pkg.
- PLTSP500Q    Spectrum IV Plate (Only) 10/Pkg.
- PLT1104      Mueller Hinton Agar; 10/Pkg. (Individually Wrapped)
- EQGLB2850   Disposable 10 ul Calibrated Inoculating Loop. 100/Pkg.
- A258062WC   Cotton Tip Sterile Applicators; 100/Pkg.
- BD4361181-5   Oxidase Mini-droppers 5/Pkg.

*\*Spectrum chromogenic agars and kits are packaged for veterinary use only. Visit [www.vetlab.com/spectrum-agar](http://www.vetlab.com/spectrum-agar) for complete specifications.*