


Detect Platelet Disorders at the Point-of-Care


Used to assess platelet function, the Buccal Mucosal Bleeding Time (BMBT) should be a part of your routine workup for any patient presented with a bleeding disorder. This is especially true for those breeds predisposed to VWD (von Willebrand's Disease). Patients with dysfunctional platelets may fail to clot even when counts are normal. The Tenderfoot® spring-loaded lancet delivers a precisely calibrated 1 mm X 2.5 mm incision allowing you to control procedural variables that could affect results.

- Tests for platelet dysfunction and adhesion.
- Aids diagnosis of VWD.
- Simple procedure applicable to any practice.
- Economical to perform.


Incise



Blot



Read Results



Complete BMBT test kit includes 10 Tenderfoot® single use automatic lancets, blotting paper, stopwatch, and instructions. Visit our website at www.vetlab.com and click on [Quicklearn](#) to view our "In-House Diagnosis of Bleeding Disorders" tutorial.

Distributed exclusively by:

Vetlab Supply · 800.330.1522 · 305.232.8421 Fax · info@vetlab.com

Tenderfoot is a registered trademark of Werfen USA