

# PowerSpin™ DX Centrifuge

UNICO® is proud to present the PowerSpin™ DX Centrifuge


With 4 rotor options this centrifuge is a great answer for every day laboratory tasks.

- 24 place fixed rotor for 2-15ml tubes (#C8724 )
- 6 place horizontal rotor for 2-15ml tubes ( #C8706 )
- 4 place horizontal rotor for up to 50ml tubes (#C8704 )
- 4 place fixed rotor for 60ml syringes (#C8760 )

The PowerSpin™ DX has variable speed 500 – 4000 RPM, to insure a clean and complete separation. Quiet and dependable, with a maintenance free motor, it can be used to process blood, urine, semen and other specimens.



Model #C8724



The all digital front panel makes it easy to select speed and run time. Oversized LCD display shows status, time, speed (RPM) and RCF (g). Precise speed control enables you to comply with CLIA recommendations for specific specimens.

- Zero RPM Safety Latching System prevents lid from being opened while rotor is spinning
- Quiet running, only 50 decibels at 3400 RPM
- 99 minute digital timer , 1 sec increment
- 500-4000 RPM, 100 RPM increment

- Remembers last setting for quick repeat runs
- Adjustable Acceleration and Deceleration
- Stainless steel rotor chamber for easy cleaning
- Suction cup rubber feet eliminate vibration
- Available in 110V and 220V - Contact UNICO for details



## UNICO Professional Digital Non-Contact Photo-Tachometer

- Highly accurate, scientific tool
- Easy to operate
- Standard adhesive reflective tape (included)
- Aim and measure--RPMs instantly display on large LCD screen.
- Includes custom protective cushioned storage case
- Order pn# C800-75