

HEATED STAGE MICROSCOPE WITH PHASE

- Complete Semen Evaluation System
- Temperature Range 25-50°C (77-122°F)
- Binocular or Video-Ready Head
- Brightfield or full Phase options
- Choice of Achromat, Semi-Plan or Plan objectives
- Halogen or LED illumination
- Digital temperature Controller
- Koehler field diaphragm for illumination control
- Digital USB or S-Video Camera options



Model G495T shown with optional digital Camera and C-Mount adapter. Contact UNICO for product details.

Model #	G390	G390SP	G390PL	G393	G393SP	G393PL	G395	G395T
Head	Binocular	Binocular	Binocular	Trinocular	Trinocular	Trinocular	Binocular	Trinocular
Objectives 4x10x40x100	Achromat	Semi-Plan	Plan	Achromat	Semi-Plan	Plan	Plan Phase 40x10x100x	Plan Phase 40x10x100x
Illumination	Halogen							
Model #	G490	G490SP	G490PL	G493	G493SP	G493PL	G495	G495T
Head	Binocular	Binocular	Binocular	Trinocular	Trinocular	Trinocular	Binocular	Trinocular
Objectives 4x10x40x100	Achromat	Semi-Plan	Plan	Achromat	Semi-Plan	Plan	Plan Phase 40x10x100x	Plan Phase 40x10x100x
Illumination	LED							