



# *MS* **Spectrum<sup>tm</sup>** **Mastitis Culture System**

## *A True Breakthrough in Mastitis Diagnosis*

Each year Mastitis results in nearly 2 billion dollars in losses to the dairy industry. In just 24 hours, the advanced Spectrum MS<sup>TM</sup> Culture System can provide a presumptive identification for many of the major organisms responsible for this economically devastating disease.



- *Formulated with Spectrum<sup>tm</sup> chromogenic agars to produce uniquely pigmented colonies.*
- *Easily differentiate infectious agents based on color and colony morphology.*
- *Includes selective agar for Staphylococcus.*
- *Individually wrapped plates for extended shelf life.*
- *Color interpretation chart provided free with first order.*

Distributed by: Vetlab Supply • Palmetto Bay, FL • 800.330.1522 • Fax 305.232.8421



[www.vetlab.com](http://www.vetlab.com)