


# Imagine the Possibilities!



The miPlatform Smart Capture System instantly transforms any smartphone into a high quality eyepiece camera. Just slide the adapter into place, align your smartphone with the adjustable grips, and start snapping away! Your captured images can be shared via SMS, email, posting to the web, or sent to a projector or computer monitor for real time viewing.



Distributed by:  
Vetlab Supply 800.330.1522 • 305.232.8421 Fax • info@vetlab.com



The revolutionary miPlatform Smart Capture System integrates the power of your existing microscope with the precision of your smartphone camera to create an essential collaborative tool with unlimited applications. Installation and alignment are simple and because there are no optical elements, the miPlatform Capture System works with any iPhone, Android, or Windows smartphone. Add captions, time/date stamp, and even print your images wirelessly with the special iPhone app available free from the Apple store.

The miPlatform universal adapter has 3 adjustable grips that lock into place to securely hold any size smartphone and an adjustable mounting collar that accommodate eyepieces up to 30 mm. The device is secured to the mounting collar with a single set screw allowing it to be easily removed for normal use and replaced in just seconds. For greater convenience, a secondary eyepiece (sold separately) can be permanently secured to the collar to create a single, easy to install unit. Placing miPlatform onto the microscope becomes just a matter of removing one of the existing eyepieces and sliding the assembly into place, virtually eliminating the need for realignment.

Visit [www.vetlab.com/miPlatform](http://www.vetlab.com/miPlatform) or call 800.330.1522 to place your order.



*Specify 23 mm or 30 mm when ordering accessory eyepieces.*

